

## **Pick up calf milk from robot (am/pm)**

**Things needed:** chemical, scrub brush.

- Go into robot room that has M4USE.
- Open the M4USE cover.
- Take the bucket out of the frame that is filled with milk.
- Bring to calf shed.
- When finished feeding calves, rinse the bucket(s) out with cold water.
- Use a drop of chemical (detergent or acid, NOT TOGETHER) and scrub the bucket(s) clean.
- Rinse the bucket(s) out with warm water.
- Return the bucket(s) back to the M4USE frame and close the cover.
- Confirm emptying the bucket for empty bucket position on robot screen by pushing "Empty".
- Close the M4USE tab.